CERTEST RHDV + RHDVb

ONE STEP *Rabbit Hemorrhagic
Disease Virus*
COMBO CARD TEST

CERTEST BIOTEC S.L.

ENGLISH**INTENDED USE**

CerTest RHDV + RHDVb one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of the "classic" and the "new variant" virions of *rabbit hemorrhagic disease virus* (RHD) in rabbit tissues samples (preferably liver) or liquid exudate from the rabbit abdominal cavity.

CerTest RHDV + RHDVb combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of rabbit hemorrhagic disease.

INTRODUCTION

Rabbit hemorrhagic disease (RHD) is an acute, highly contagious and fatal disease of the European rabbit (*Oryctolagus cuniculus*), caused by a calicivirus (genus *Lagovirus*).

Infection can occur by nasal, conjunctival or orally. RHD transmission is facilitated by the high stability of the virus in the environment. The incubation period ranges from 1 to 3 days, and the death usually occurs 12–36 hours after the onset of fever. Clinical signs have been described mainly in the acute infection (nervous and respiratory signs, apathy and anorexia). Specific and clear lesions appear both macroscopic and microscopic. There is a primary hepatic necrosis and a massive disseminated intravascular coagulopathy in all organs and tissues. The most severe lesions are found in the liver, trachea and lungs. In almost every organs petechiae are demonstrated, which are accompanied by poor blood coagulation. In rabbits under 4-5 weeks of age, infection with the "classic" RHDV is subclinical.

Recently a "new variant" virus (RHDVb) which affects young rabbits and farmed adult animals previously immunized with vaccines against "classic" RHDV has been detected.

TEST PRINCIPLE

CerTest RHDV + RHDVb is based on the principle of a qualitative immunochromatographic assay for the determination of the "classic" and the "new variant" virions of *rabbit hemorrhagic disease virus* in rabbit tissues samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against virus like particles (VLP's) and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-RHDV) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against virus like particles (VLP's) and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-RHDVb) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is RHDV positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-RHDV monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is RHDVb positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-RHDVb monoclonal antibodies-red polystyrene microspheres) in strip B which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-RHDV antibodies present on the membrane of strip A (test line) and the anti-RHDVb antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no RHDV presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-RHDV and the anti-RHDVb antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest RHDV + RHDVb combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Tissue samples (preferably liver) or abdominal liquid exudate should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) prior to testing. For longer storage, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing.


Tissue specimen preparation (see illustration 1):

1. Take approximately 0.1g of liver tissue, mix with 1ml PBS (or tap water) and homogenize with the aid of a disposable pipette. Leave the homogenate to stand to allow the particulate material to settle and use the supernatant (approx. 10% tissue homogenate) in further steps (1).
2. Take out the cap of the collection tube (2)
3. Add 7 drops of the tissue homogenate supernatant using a dispensable pipette, in the collection tube (3).
4. Close the tube with the diluent and sample. Shake the tube in order to assure good sample dispersion (4).

Abdominal liquid exudate specimen preparation (see illustration 1):

1. Take out the cap of the collection tube (2)
2. Abdominal liquid exudate samples can be used directly without previous treatment. Add 7 drops of the abdominal liquid exudate using a dispensable pipette, in the collection tube (3).
3. Close the tube with the diluent and sample. Shake the tube in order to assure good sample dispersion (4).

Illustration 1

MATERIALS

MATERIALS PROVIDED

- CerTest RHDV + RHDVb combo card tests
- Instructions for use
- Disposable pipettes
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE



Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the collection tube in order to assure good sample dispersion.
2. Remove CerTest RHDV + RHDVb combo card test from its sealed bag just before using it.
3. Take the collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter A (5) and 4 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.

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INTERPRETATION OF THE RESULTS (please refer to the illustration below)

- A: Green→Negative RHDV
B: Green→Negative RHDVb
- A and B: Green/Red→Positive RHDVb
- A: Green/Red→Positive RHDV
B: Green →Negative RHDVb

INVALID RESULTS

- A: Green→Negative RHDV
B: Green/Red→Positive RHDVb
-

	A (RHDV)	B (RHDVb)	Interpretation of the results
1.	- GREEN	- GREEN	There are no "classic" (RHDV) and "new variant" (RHDVb) virions of <i>rabbit hemorrhagic disease virus</i> .
2.	+ GREEN-RED	+ GREEN-RED	There is presence of virions of "new variant" (RHDVb) of <i>rabbit hemorrhagic disease virus</i> . To dismiss a very unusual case of co-infection with other type of RHDV, some complementary RT-PCR assays should be made.
3.	+ GREEN-RED	- GREEN	There is "classic" (RHDV) virions of <i>rabbit hemorrhagic disease virus</i> .
4.	- GREEN	+ GREEN-RED	Invalid result: If B is positive, A has to be also positive, we recommend repeating the assay using the same sample with another test. If the same result is obtained, some complementary RT-PCR assays should be made to verify the possible type of RHDV present.
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

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NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of virions present in the specimen. However, neither the quantitative value nor the rate of increase in virions can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high virus concentration to faint when the virus concentration is close to the detection limit value of the test.
4. CerTest RHDV + RHDVb should be used from tissue or exudate samples from rabbits. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper specimens must be obtained.
5. Positive results determine the presence of the "classic" and/or the "new variant" virions of isolate of *rabbit hemorrhagic disease virus* in tissue sample or exudate. A positive result should be followed up with other diagnostics to confirm the results. A confirmed infection should only be made by a specialised professional after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A very low number of negative results have been observed with the CerTest RHDV + RHDVb test (strip A RHDV) when the responsible virus belongs to the Genogroup 6 (RHDVa), although other similar isolates (Hartmannsdorf, Triptis or Rossi) are detected efficiently with this test. In the case that suspected *rabbit haemorrhagic disease virus* clinical and laboratory findings exist and a negative result is obtained with the CerTest RHDV+RHDVb test, an ELISA + RT-PCR assay should be carried out to discard the presence of RHDVa (genogroup 6).
7. A negative result is not meaningful because of it is possible the amount of virions in the sample is lower than the detection limit value. In the case that suspected *rabbit haemorrhagic disease virus* clinical and laboratory findings exist and a negative result is obtained with the CerTest RHDV test, additional assays (ELISA + RT-PCR) should be carried out to confirm the result.

EXPECTED VALUES

Rabbit Hemorrhagic Disease is characterized by high morbidity and high mortality (70-90%), and spreads very rapidly through direct and indirect transmission.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity (RHDV + RHDVb)

An evaluation, with positive and negative samples, was performed comparing the results obtained by an immunochromatographic test (CerTest RHDV + RHDVb, CerTest) with positive and negative rabbit liver specimens; these samples were evaluated compared with an ELISA assay (INGEZIM RHDV DAS (Ingenasa)). The results were as follows:

		ELISA + RT-PCR (RHDV)		
		+	-	Total
IC test: CerTest RHDV + RHDVb (RHDV)	+	69	0	69
	-	2 (*)	10	12
	Total	71	10	81

	Sensitivity	Specificity	PPV	NPV
CerTest RHDV + RHDVb vs ELISA + RT-PCR (RHDV)	>97%	>99%	>99%	>83%

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(*) The two positive samples which, being positive in the ELISA assay + RT-PCR (RHDV), were negative in the immunochromatographic test CerTest RHDV (**Strip A**) belong to the genogroup 6 and are an antigenic subtype known as RHDVa.

		ELISA + RT-PCR (RHDVb)		
		+	-	Total
IC test: CerTest RHDV + RHDVb (RHDVb)	+	62	0	62
	-	0	19	19
	Total	62	19	81

	Sensitivity	Specificity	PPV	NPV
CerTest RHDV + RHDVb vs ELISA + RT-PCR (RHDVb)	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect RHDV (RHDV classical and RHDVb new variant) of rabbit hemorrhagic disease virus using CerTest RHDV + RHDVb.

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